



PKHG1 Polyclonal Antibody

Catalog No	YP-Ab-05938
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	PLEKHG1 KIAA1209
Immunogen	Synthesized peptide derived from human protein . at AA range: 780-860
Specificity	PKHG1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PKHG1;PLEKHG1 KIAA1209;Pleckstrin homology domain-containing family G member 1
Observed Band	152kD
Calculated Molecular Weight	152kD
Cell Pathway	
Tissue Specificity	Brain,Epithelium,Lymph,Pooled,Testicle,
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Contains 1 DH (DBL-homology) domain.,similarity:Contains 1 PH domain.,
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Contains 1 DH (DBL-homology) domain.,similarity:Contains 1 PH domain.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images