



PCF11 Polyclonal Antibody

Catalog No	YP-Ab-05967
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	PCF11 KIAA0824
Immunogen	Synthesized peptide derived from human protein . at AA range: 490-570
Specificity	PCF11 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA0824;Pre-mRNA cleavage complex 2 protein Pcf11;Pre-mRNA cleavage complex II protein Pcf11;PCF11;PCF11 KIAA0824;Pre-mRNA cleavage complex 2 protein Pcf11 (Pre-mRNA cleavage complex II protein Pcf11)
Observed Band	173 kDa
Calculated Molecular Weight	1555 aa, 173 kDa
Cell Pathway	Nucleus.
Tissue Specificity	Brain,Epithelium,Testis,Uterus,
Function	function:Component of pre-mRNA cleavage complex II.,sequence caution:Contaminating sequence. Potential poly-A sequence.,similarity:Contains 1 CID domain.,
Background	The protein encoded by this gene binds to CLP1 to form pre-mRNA cleavage factor IIm. The encoded protein is necessary for efficient Pol II transcription termination and may be involved in degradation of the 3' product of polyA site cleavage. [provided by RefSeq, Oct 2016],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images