



CLP1 Polyclonal Antibody

Catalog No	YP-Ab-05968
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	CLP1 HEAB
Immunogen	Synthesized peptide derived from human protein . at AA range: 50-130
Specificity	CLP1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CLP1;hClp1;HEAB;Polynucleotide kinase Clp1;CLP1 HEAB;Polyribonucleotide 5'-hydroxyl-kinase Clp1 (EC 2.7.1.78) (Polynucleotide kinase Clp1) (Pre-mRNA cleavage complex II protein Clp1)
Observed Band	46kD
Calculated Molecular Weight	46kD
Cell Pathway	Nucleus .
Tissue Specificity	Lung,Placenta,
Function	catalytic activity:ATP + 5'-dephospho-DNA = ADP + 5'-phospho-DNA.,catalytic activity:ATP + 5'-dephospho-RNA = ADP + 5'-phospho-RNA., cofactor:Magnesium, manganese or nickel.,function:Polynucleotide kinase that can phosphorylate the 5'-hydroxyl groups of double-stranded RNA (dsRNA), single-stranded RNA (ssRNA), double stranded DNA (dsDNA) and double-stranded DNA:RNA hybrids. dsRNA is phosphorylated more efficiently than dsDNA, and the RNA component of a DNA:RNA hybrid is phosphorylated more efficiently than the DNA component. Appears to have roles in both tRNA splicing and mRNA 3'-end formation. Component of the tRNA splicing endonuclease complex. Phosphorylates the 5'-terminus of the tRNA 3'-exon during tRNA splicing; this phosphorylation event is a prerequisite for the subsequent ligation



of the two exon halves and the production of a mature tRNA. Component of the pre-mRNA cleavage compl

Background

This gene encodes a member of the Clp1 family. The encoded protein is a multifunctional kinase which is a component of the tRNA splicing endonuclease complex and a component of the pre-mRNA cleavage complex II. This protein is implicated in tRNA, mRNA, and siRNA maturation. Mutations in this gene are associated with pontocerebellar hypoplasia type 10 (PCH10). Alternatively splice transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2014],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night