



POTEE Polyclonal Antibody

Catalog No	YP-Ab-05978
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	POTEE A26C1A POTE2
Immunogen	Synthesized peptide derived from human protein . at AA range: 350-430
Specificity	POTEE Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	POTEE;POTEE A26C1A POTE2;POTE ankyrin domain family member E (ANKRD26-like family C member 1A) (Prostate;ovary;testis-expressed protein on chromosome 2) (POTE-2)
Observed Band	118kD
Calculated Molecular Weight	118kD
Cell Pathway	extracellular space,extracellular exosome,blood microparticle,
Tissue Specificity	Hepatocellular carcinoma,Prostate,
Function	similarity:Belongs to the actin family.,similarity:Contains 5 ANK repeats.,similarity:In the C-terminal section; belongs to the actin family.,similarity:In the N-terminal section; belongs to the POTE family.,
Background	similarity:Belongs to the actin family.,similarity:Contains 5 ANK repeats.,similarity:In the C-terminal section; belongs to the actin family.,similarity:In the N-terminal section; belongs to the POTE family.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images