



FOH1B Polyclonal Antibody

Catalog No	YP-Ab-05979
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	FOLH1B PSMAL GIG26
Immunogen	Synthesized peptide derived from human protein . at AA range: 370-450
Specificity	FOH1B Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FOH1B;FOLH1B PSMAL GIG26;Putative N-acetylated-alpha-linked acidic dipeptidase (NAALADase) (EC 3.4.-.-) (Cell growth-inhibiting gene 26 protein) (Prostate-specific membrane antigen-like protein) (Putative folate hydrolase 1B)
Observed Band	48kD
Calculated Molecular Weight	48kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Kidney and liver. Not expressed in the prostate.
Function	proteolysis,
Background	FOLH1B (Folate Hydrolase 1B) is a Protein Coding gene. Among its related pathways are Metabolism and Amino acid synthesis and interconversion (transamination). Has both folate hydrolase and N-acetylated-alpha-linked-acidic dipeptidase (NAALADase) activity. Exhibits a dipeptidyl-peptidase IV type activity.
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images