



PPM1H Polyclonal Antibody

Catalog No	YP-Ab-05984
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	PPM1H ARHCL1 KIAA1157 URCC2
Immunogen	Synthesized peptide derived from human protein . at AA range: 290-370
Specificity	PPM1H Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ARHCL1;EC:3.1.3.16;KIAA1157;Protein phosphatase 1H;URCC2;PPM1H;PPM1H ARHCL1 KIAA1157 URCC2;Protein phosphatase 1H (EC 3.1.3.16)
Observed Band	50-56 kDa
Calculated Molecular Weight	82 aa, 9 kDa
Cell Pathway	Nucleus . Cytoplasm .
Tissue Specificity	Brain,
Function	catalytic activity:A phosphoprotein + H(2)O = a protein + phosphate., similarity:Belongs to the PP2C family.,similarity:Contains 1 PP2C-like domain.,
Background	catalytic activity:A phosphoprotein + H(2)O = a protein + phosphate., similarity:Belongs to the PP2C family.,similarity:Contains 1 PP2C-like domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images