



S1A7A Polyclonal Antibody

Catalog No	YP-Ab-05992
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	S100A7A S100A15 S100A7L1
Immunogen	Synthesized peptide derived from human protein . at AA range: 30-110
Specificity	S1A7A Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	S1A7A;S100A7A S100A15 S100A7L1;Protein S100-A7A (S100 calcium-binding protein A15) (S100 calcium-binding protein A7-like 1) (S100 calcium-binding protein A7A)
Observed Band	11kD
Calculated Molecular Weight	11kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Overexpressed in psoriasis.
Function	function:May be involved in epidermal differentiation and inflammation and might therefore be important for the pathogenesis of psoriasis and other diseases., similarity:Belongs to the S-100 family.,similarity:Contains 2 EF-hand domains., tissue specificity:Overexpressed in psoriasis.,
Background	function:May be involved in epidermal differentiation and inflammation and might therefore be important for the pathogenesis of psoriasis and other diseases., similarity:Belongs to the S-100 family.,similarity:Contains 2 EF-hand domains., tissue specificity:Overexpressed in psoriasis.,

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images