



DEN5A Polyclonal Antibody

Catalog No	YP-Ab-06005
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	DENND5A KIAA1091 RAB6IP1
Immunogen	Synthesized peptide derived from human protein . at AA range: 650-730
Specificity	DEN5A Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DEN5A_HUMAN;DENN domain-containing protein 5A;DENN/MADD domain containing 5A;Dennd5a;FLJ22354;FLJ33829;FLJ43455;KIAA191;RAB6 interacting protein 1;Rab6-interacting protein 1;Rab6IP1;DEN5A;DENND5A KIAA1091 RAB6IP1;DENN domain-containing protein 5A (Rab6-interacting protein 1) (Rab6IP1)
Observed Band	141kD
Calculated Molecular Weight	147kD(Calculated).
Cell Pathway	Golgi apparatus membrane .
Tissue Specificity	Brain,Placenta,Skin,Uterus,
Function	function:May be involved in Rab6-mediated GTPase signaling., similarity:Belongs to the RAB6IP1 family.,similarity:Contains 1 dDENN domain., similarity:Contains 1 DENN domain.,similarity:Contains 1 PLAT domain., similarity:Contains 1 uDENN domain.,similarity:Contains 2 RUN domains., subunit:Probably interacts with RAB6A bound to GTP.,
Background	function:May be involved in Rab6-mediated GTPase signaling., similarity:Belongs to the RAB6IP1 family.,similarity:Contains 1 dDENN domain., similarity:Contains 1 DENN domain.,similarity:Contains 1 PLAT domain.,

similarity:Contains 1 uDENN domain.,similarity:Contains 2 RUN domains.,
subunit:Probably interacts with RAB6A bound to GTP.,

**matters needing
attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images