



RGL1 Polyclonal Antibody

Catalog No	YP-Ab-06013
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	RGL1 KIAA0959 RGL
Immunogen	Synthesized peptide derived from human protein . at AA range: 210-290
Specificity	RGL1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Rho-type GTPase-activating protein 9;Rho GTPase-activating protein 9;RGL1; 10C;RGL1 KIAA0959 RGL;Ral guanine nucleotide dissociation stimulator-like 1 (RalGDS-like 1)
Observed Band	70 kDa, 65 kDa
Calculated Molecular Weight	750 aa, 83 kDa
Cell Pathway	intracellular,cytosol,
Tissue Specificity	Expressed in a wide variety of tissues with strong expression being seen in the heart, brain, kidney, spleen and testis.
Function	function:Probable guanine nucleotide exchange factor.,similarity:Contains 1 N-terminal Ras-GEF domain.,similarity:Contains 1 Ras-associating domain.,similarity:Contains 1 Ras-GEF domain.,subunit:Interacts with Ras.,tissue specificity:Expressed in a wide variety of tissues with strong expression being seen in the heart, brain, kidney, spleen and testis.,
Background	function:Probable guanine nucleotide exchange factor.,similarity:Contains 1 N-terminal Ras-GEF domain.,similarity:Contains 1 Ras-associating domain.,similarity:Contains 1 Ras-GEF domain.,subunit:Interacts with Ras.,tissue specificity:Expressed in a wide variety of tissues with strong expression being

seen in the heart, brain, kidney, spleen and testis.,

**matters needing
attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images