




RAB30 Polyclonal Antibody

Catalog No	YP-Ab-06046
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	RAB30
Immunogen	Synthesized peptide derived from human protein . at AA range: 80-160
Specificity	RAB30 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RAB30;Ras related protein Rab 30;Ras-related protein Rab-30
Observed Band	22kD
Calculated Molecular Weight	22kD
Cell Pathway	Membrane ; Lipid-anchor ; Cytoplasmic side . Golgi apparatus, trans-Golgi network . Cytoplasm . Golgi apparatus .
Tissue Specificity	B-cell,Brain,Cerebellum,Human uterus,Melanocyte,
Function	similarity:Belongs to the small GTPase superfamily. Rab family.,
Background	similarity:Belongs to the small GTPase superfamily. Rab family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night