




RAB36 Polyclonal Antibody

Catalog No	YP-Ab-06049
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	RAB36
Immunogen	Synthesized peptide derived from human protein . at AA range: 170-250
Specificity	RAB36 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RAB 36;Rab36;RAB36 member RAS oncogene family;RAB36_HUMAN;Ras related protein Rab 36;Ras related protein Rab36;Ras-related protein Rab-36; Small GTP binding protein Rab36
Observed Band	36kD
Calculated Molecular Weight	36kD
Cell Pathway	Golgi apparatus membrane; Lipid-anchor.
Tissue Specificity	Ubiquitously present in all tissues examined.
Function	function:Protein transport. Probably involved in vesicular traffic., similarity:Belongs to the small GTPase superfamily. Rab family.,tissue specificity:Ubiquitously present in all tissues examined.,
Background	function:Protein transport. Probably involved in vesicular traffic., similarity:Belongs to the small GTPase superfamily. Rab family.,tissue specificity:Ubiquitously present in all tissues examined.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night