



RAB6B Polyclonal Antibody

Catalog No	YP-Ab-06057
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	RAB6B
Immunogen	Synthesized peptide derived from human protein . at AA range: 80-160
Specificity	RAB6B Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:3.6.5.2;Ras related protein Rab 6B;Ras-related protein Rab-6B;RAB6B
Observed Band	22kD
Calculated Molecular Weight	22kD
Cell Pathway	Golgi apparatus membrane ; Lipid-anchor . Endoplasmic reticulum-Golgi intermediate compartment . Cytoplasmic vesicle . Colocalizes with BICD1 at vesicular structures that align along microtubules. .
Tissue Specificity	Predominantly expressed in brain.
Function	function:Seems to have a role in retrograde membrane traffic at the level of the Golgi complex.,similarity:Belongs to the small GTPase superfamily. Rab family., subunit:Interacts with RAB6KIFL.,tissue specificity:Predominantly expressed in brain.,
Background	function:Seems to have a role in retrograde membrane traffic at the level of the Golgi complex.,similarity:Belongs to the small GTPase superfamily. Rab family., subunit:Interacts with RAB6KIFL.,tissue specificity:Predominantly expressed in brain.,


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from A431 cells, primary antibody was diluted at 1:1000, 4° over night