



RHG30 Polyclonal Antibody

Catalog No	YP-Ab-06088
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	ARHGAP30
Immunogen	Synthesized peptide derived from human protein . at AA range: 230-310
Specificity	RHG30 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RHG30;ARHGAP30;Rho GTPase-activating protein 30 (Rho-type GTPase-activating protein 30)
Observed Band	121kD
Calculated Molecular Weight	121kD
Cell Pathway	Cytoplasmic vesicle .
Tissue Specificity	Blood,Endometrium,Epithelium,Spleen,T-cell,Testis,Thymus,
Function	function:GTPase activator for the Rho-type GTPases by converting them to an inactive GDP-bound state.,similarity:Contains 1 Rho-GAP domain.,
Background	function:GTPase activator for the Rho-type GTPases by converting them to an inactive GDP-bound state.,similarity:Contains 1 Rho-GAP domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
