



RL1D1 Polyclonal Antibody

Catalog No	YP-Ab-06105
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	RSL1D1 CATX11 CSIG PBK1 L12
Immunogen	Synthesized peptide derived from human protein . at AA range: 280-360
Specificity	RL1D1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CATX 11;CATX11;CSIG;L12;PBK1;Protein PBK1;RSL1D1;RL1D1;RSL1D1 CATX11 CSIG PBK1 L12;Ribosomal L1 domain-containing protein 1 (CATX-11) (Cellular senescence-inhibited gene protein) (Protein PBK1)
Observed Band	55-70 kDa
Calculated Molecular Weight	490 aa, 55 kDa
Cell Pathway	Nucleus, nucleolus . Colocalizes with ING1 in the nucleolus after UV stress. .
Tissue Specificity	Expressed at high intensities in the heart, skeletal muscle, and placenta.
Function	similarity:Belongs to the ribosomal protein L1P family. Highly divergent.,tissue specificity:Placenta.,
Background	similarity:Belongs to the ribosomal protein L1P family. Highly divergent.,tissue specificity:Placenta.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
