



BMS1 Polyclonal Antibody

Catalog No	YP-Ab-06107
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	BMS1 BMS1L KIAA0187
Immunogen	Synthesized peptide derived from human protein . at AA range: 210-290
Specificity	BMS1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BMS1;BMS1 BMS1L KIAA0187;Ribosome biogenesis protein BMS1 homolog (Ribosome assembly protein BMS1 homolog)
Observed Band	141kD
Calculated Molecular Weight	141kD
Cell Pathway	Nucleus, nucleolus .
Tissue Specificity	Bone marrow,Clones donated by Kazusa DNA Research Inst.,Skin,
Function	function:May act as a molecular switch during maturation of the 40S ribosomal subunit in the nucleolus.,similarity:Belongs to the BMS1/TSR1 family. BMS1 subfamily.,
Background	This gene likely encodes a ribosome assembly protein. A similar protein in yeast functions in 35S-rRNA processing, which includes a series of cleavage steps critical for formation of 40S ribosomes. Related pseudogenes exist on chromosomes 2, 9, 10, 15, 16, and 22.[provided by RefSeq, Mar 2009],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images