



RPF1 Polyclonal Antibody


Catalog No	YP-Ab-06109
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	RPF1 BXDC5
Immunogen	Synthesized peptide derived from human protein . at AA range: 180-260
Specificity	RPF1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Brix domain-containing protein 5;Ribosome biogenesis protein RPF1;Ribosome production factor 1;Ribosome Production Factor 1 Homolog;RPF1;RPF1 BXDC5; Ribosome production factor 1 (Brix domain-containing protein 5) (Ribosome biogenesis protein RPF1)
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	Nucleus, nucleolus .
Tissue Specificity	Amygdala,B-cell,Teratocarcinoma,
Function	function:May be required for ribosome biogenesis.,similarity:Contains 1 Brix domain.,
Background	function:May be required for ribosome biogenesis.,similarity:Contains 1 Brix domain.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from DU145 cells,
primary antibody was diluted at 1:1000, 4° over night