



RPAP3 Polyclonal Antibody

Catalog No	YP-Ab-06114
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	RPAP3
Immunogen	Synthesized peptide derived from human protein . at AA range: 280-360
Specificity	RPAP3 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RNA polymerase II-associated protein 3;RPAP3
Observed Band	76 kDa
Calculated Molecular Weight	76 kDa
Cell Pathway	R2TP complex,
Tissue Specificity	Blood,Epithelium,Tongue,
Function	function:Forms an interface between the RNA polymerase II enzyme and chaperone/scaffolding protein, suggesting that it is required to connect RNA polymerase II to regulators of protein complex formation.,similarity:Belongs to the RPAP3 family.,similarity:Contains 7 TPR repeats.,subunit:Tightly associated with the RNA polymerase II complex.,
Background	This gene encodes an RNA polymerase II-associated protein. The encoded protein may function in transcriptional regulation and may also regulate apoptosis. Alternatively spliced transcript variants have been described. [provided by RefSeq, Mar 2009],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images