




PRS23 Polyclonal Antibody

Catalog No	YP-Ab-06123
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	PRSS23 ZSIG13 UNQ270/PRO307
Immunogen	Synthesized peptide derived from human protein . at AA range: 150-230
Specificity	PRS23 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PRS23;PRSS23 ZSIG13 UNQ270/PRO307;Serine protease 23 (EC 3.4.21.-) (Putative secreted protein Zsig13)
Observed Band	42kD
Calculated Molecular Weight	42kD
Cell Pathway	Secreted .
Tissue Specificity	Cervix,Umbilical vein,Uterus,
Function	similarity:Belongs to the peptidase S1 family.,
Background	This gene encodes a conserved member of the trypsin family of serine proteases. Mouse studies found a decrease of mRNA levels of this gene after ovulation was induced. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2014],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A431 cells,
primary antibody was diluted at 1:1000, 4°over night