



PP2BC Polyclonal Antibody

Catalog No	YP-Ab-06152
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	PPP3CC CALNA3 CNA3
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	PP2BC Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CNA3;CAM-PRP catalytic subunit;CAM PRP catalytic subunit;CALNA3; Calmodulin-dependent calcineurin A subunit gamma isoform;PP2BC;PPP3CC CALNA3 CNA3;Serine/threonine-protein phosphatase 2B catalytic subunit gamma isoform (EC 3.1.3.16) (CAM-PRP catalytic subunit) (Calcineurin;testis-specific catalytic subunit) (Calmodulin-dependent calcineurin A s
Observed Band	57 kDa
Calculated Molecular Weight	57 kDa
Cell Pathway	Mitochondrion . Localizes in the mitochondria in a SPATA33-dependent manner.
Tissue Specificity	Testis.
Function	catalytic activity:A phosphoprotein + H(2)O = a protein + phosphate., cofactor: Binds 1 Fe(3+) ion per subunit.,cofactor: Binds 1 zinc ion per subunit., function:Calcium-dependent, calmodulin-stimulated protein phosphatase. This subunit may have a role in the calmodulin activation of calcineurin., similarity:Belongs to the PPP phosphatase family. PP-2B subfamily., subunit:Composed of two components (A and B), the A component is the catalytic subunit and the B component confers calcium sensitivity.,tissue specificity:Testis.,

Background

Calcineurin is a calcium-dependent, calmodulin-stimulated protein phosphatase involved in the downstream regulation of dopaminergic signal transduction. Calcineurin is composed of a regulatory subunit and a catalytic subunit. The protein encoded by this gene represents one of the regulatory subunits that has been found for calcineurin. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2011],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images