



SCAFB Polyclonal Antibody

Catalog No	YP-Ab-06161
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	SCAF11 CASP11 SFRS2IP SIP1 SRSF2IP
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	SCAFB Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CASP11;CTD-associated SR protein 11;NY-REN-40;Protein SCAF11;Renal carcinoma antigen NY-REN-40;SCAFB;SCAF11 CASP11 SFRS2IP SIP1 SRSF2IP;Protein SCAF11 (CTD-associated SR protein 11) (Renal carcinoma antigen NY-REN-40) (SC35-interacting protein 1) (SR-related and CTD-associated factor 11) (SRSF2-interacting protein) (Serine/arginine-ri
Observed Band	240 kDa
Calculated Molecular Weight	1463 aa, 165 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Widely expressed.
Function	function:Plays a role in pre-mRNA alternative splicing by regulating spliceosome assembly.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., sequence caution:Contaminating sequence. Potential poly-A sequence., subunit:Interacts with SFRS2, U2AF2 and SNRNP70.,tissue specificity:Widely expressed.,
Background	function:Plays a role in pre-mRNA alternative splicing by regulating spliceosome assembly.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., sequence caution:Contaminating sequence. Potential poly-A sequence.,

subunit:Interacts with SFRS2, U2AF2 and SNRNP70.,tissue specificity:Widely expressed.,

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images