



SH23A Polyclonal Antibody

Catalog No	YP-Ab-06164
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	SH2D3A NSP1 UNQ175/PRO201
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	SH23A Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	UNQ175/PRO201;SH2 domain-containing protein 3A;NSP1;Novel SH2-containing protein 1;SH23A;SH2D3A NSP1 UNQ175/PRO201;SH2 domain-containing protein 3A (Novel SH2-containing protein 1)
Observed Band	64 kDa
Calculated Molecular Weight	576 aa, 63 kDa
Cell Pathway	intracellular,
Tissue Specificity	Weakly expressed in placenta, fetal kidney, fetal lung, adult pancreas, adult kidney and adult lung.
Function	function:May play a role in JNK activation.,PTM:Phosphorylated on tyrosine.,similarity:Contains 1 SH2 domain.,subunit:Interacts with BCAR1.,tissue specificity:Weakly expressed in placenta, fetal kidney, fetal lung, adult pancreas, adult kidney and adult lung.,
Background	function:May play a role in JNK activation.,PTM:Phosphorylated on tyrosine.,similarity:Contains 1 SH2 domain.,subunit:Interacts with BCAR1.,tissue specificity:Weakly expressed in placenta, fetal kidney, fetal lung, adult pancreas, adult kidney and adult lung.,

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images