




SASH3 Polyclonal Antibody

Catalog No	YP-Ab-06178
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	SASH3 CXorf9 SLY
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	SASH3 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CXorf9;SAM and SH3 domain-containing protein 3;SH3 protein expressed in lymphocytes homolog;SLY;SASH3;SASH3 CXorf9 SLY;SAM and SH3 domain-containing protein 3 (SH3 protein expressed in lymphocytes homolog)
Observed Band	41kD
Calculated Molecular Weight	41kD
Cell Pathway	nucleus,cytoplasm,
Tissue Specificity	Spleen,Testis,
Function	function:May function as a signaling adapter protein in lymphocytes., similarity:Contains 1 SAM (sterile alpha motif) domain.,similarity:Contains 1 SH3 domain.,
Background	The protein encoded by this gene contains a Src homology-3 (SH3) domain and a sterile alpha motif (SAM), both of which are found in proteins involved in cell signaling. This protein may function as a signaling adapter protein in lymphocytes.[provided by RefSeq, Sep 2009],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4° over night