



NAGT1 Polyclonal Antibody

Catalog No	YP-Ab-06216
Isotype	IgG
Reactivity	Human;Rat
Applications	WB;ELISA
Gene Name	NAGLT1 KIAA1919 HSPC100
Immunogen	Synthesized peptide derived from human protein . at AA range: 400-480
Specificity	NAGT1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NAGT1;NAGLT1 KIAA1919 HSPC100;Sodium-dependent glucose transporter 1
Observed Band	56kD
Calculated Molecular Weight	56kD
Cell Pathway	Apical cell membrane ; Multi-pass membrane protein .
Tissue Specificity	Brain,Testis,Umbilical cord blood,
Function	function:May function as a sodium-dependent glucose transporter., similarity:Belongs to the major facilitator superfamily.,
Background	function:May function as a sodium-dependent glucose transporter., similarity:Belongs to the major facilitator superfamily.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

