



S20A1 Polyclonal Antibody

Catalog No	YP-Ab-06220
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	SLC20A1 GLVR1 PIT1
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	S20A1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Gibbon ape leukemia virus receptor 1;GLVR1;GLVR-1;Leukemia virus receptor 1 homolog;Phosphate transporter 1;S20A1;SLC20A1 GLVR1 PIT1;Sodium-dependent phosphate transporter 1 (Gibbon ape leukemia virus receptor 1) (GLVR-1) (Leukemia virus receptor 1 homolog) (Phosphate transporter 1) (PiT-1) (Solute carrier family 20 member 1)
Observed Band	85 kDa
Calculated Molecular Weight	679 aa, 74 kDa
Cell Pathway	Membrane; Multi-pass membrane protein.
Tissue Specificity	Ubiquitously expressed.
Function	domain:Region A confers human cells susceptibility to infection by Gibbon Ape Leukemia Virus (GaLV) and Feline leukemia virus subgroup B (FeLV-B). Substitution of Human SLC20A1 region A by region A of murine SLC20A1 prevents viral infection.,function:Sodium-phosphate symporter which plays a fundamental housekeeping role in phosphate transport, such as absorbing phosphate from interstitial fluid for normal cellular functions such as cellular metabolism, signal transduction, and nucleic acid and lipid synthesis. May play a role in extracellular matrix and cartilage calcification as well as in vascular calcification. May function as a retroviral receptor as it confers human cells



susceptibility to infection to Gibbon Ape Leukemia Virus (GaLV), Simian sarcoma-associated virus (SSAV) and Feline leukemia virus subgroup B (FeLV-B) as well as 10A1 murine leukemia virus (10A1 MLV).,induction:By p

Background

The protein encoded by this gene is a sodium-phosphate symporter that absorbs phosphate from interstitial fluid for use in cellular functions such as metabolism, signal transduction, and nucleic acid and lipid synthesis. The encoded protein is also a retroviral receptor, causing human cells to be susceptible to infection by gibbon ape leukemia virus, simian sarcoma-associated virus, feline leukemia virus subgroup B, and 10A1 murine leukemia virus.[provided by RefSeq, Mar 2011],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images