



SRSF9 Polyclonal Antibody

Catalog No	YP-Ab-06236
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	SRSF9 SFRS9 SRP30C
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	SRSF9 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SRSF9;Pre-mRNA-splicing factor SRp30C;Serine/arginine-rich splicing factor 9; Splicing factor;arginine/serine-rich 9;SRP30C;SRSF9 SFRS9 SRP30C;Serine/arginine-rich splicing factor 9 (Pre-mRNA-splicing factor SRp30C) (Splicing factor;arginine/serine-rich 9)
Observed Band	24kD
Calculated Molecular Weight	24kD
Cell Pathway	Nucleus . Cellular stresses such as heat shock may induce localization to discrete nuclear bodies termed SAM68 nuclear bodies (SNBs), HAP bodies, or stress bodies. Numerous splicing factors including SRSF1/SFRS1/SF2, SRSF7/SFRS7, SAFB and KHDRBS1/SAM68 accumulate at these structures, which may participate in the post-transcriptional regulation of mRNAs in stressed cells.
Tissue Specificity	Expressed at high levels in the heart, kidney, pancreas and placenta, and at lower levels in the brain, liver, lung and skeletal muscle.
Function	function:Plays a role in constitutive splicing and can modulate the selection of alternative splice sites.,PTM:Extensively phosphorylated on serine residues in the RS domain.,similarity:Belongs to the splicing factor SR family.,similarity:Contains 2 RRM (RNA recognition motif) domains.,subcellular location:Cellular stresses such as heat shock may induce localization to discrete



nuclear bodies termed SAM68 nuclear bodies (SNBs), HAP bodies, or stress bodies. Numerous splicing factors including SFRS1/SF2/ASF, SFRS7/9G8, SAFB and KHDRBS1/SAM68 accumulate at these structures, which may participate in the post-transcriptional regulation of mRNAs in stressed cells., subunit: Interacts with KHDRBS3 (By similarity). Interacts with NOL3/ARC/ NOP30, NSEP1/YB-1/YB1, SAFB/SAFB1, SFRS6/SRP55 and TRA2B/SFRS10. May also interact with DUSP11/PIR1., tissue specificity: Expressed at high levels in the heart, k

Background

The protein encoded by this gene is a member of the serine/arginine (SR)-rich family of pre-mRNA splicing factors, which constitute part of the spliceosome. Each of these factors contains an RNA recognition motif (RRM) for binding RNA and an RS domain for binding other proteins. The RS domain is rich in serine and arginine residues and facilitates interaction between different SR splicing factors. In addition to being critical for mRNA splicing, the SR proteins have also been shown to be involved in mRNA export from the nucleus and in translation. Two pseudogenes, one on chromosome 15 and the other on chromosome 21, have been found for this gene. [provided by RefSeq, Sep 2010],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night