



SYT14 Polyclonal Antibody

Catalog No	YP-Ab-06256
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	SYT14
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	SYT14 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Synaptotagmin XIV;SYT14;SytXIV;Synaptotagmin-14 (Synaptotagmin XIV) (SytXIV)
Observed Band	62-65 kDa
Calculated Molecular Weight	62 kDa
Cell Pathway	Membrane ; Single-pass type III membrane protein . Localized in perinuclear and submembranous regions.
Tissue Specificity	Highly expressed in fetal and adult brain tissue.
Function	function:May be involved in the trafficking and exocytosis of secretory vesicles in non-neuronal tissues. Is Ca(2+)-independent.,similarity:Belongs to the synaptotagmin family.,similarity:Contains 2 C2 domains.,subunit:Homodimer. Can also form heterodimers.,
Background	synaptotagmin 14(SYT14) Homo sapiens This gene is a member of the synaptotagmin gene family and encodes a protein similar to other family members that mediate membrane trafficking in synaptic transmission. The encoded protein is a calcium-independent synaptotagmin. Mutations in this gene are a cause of autosomal recessive spinocerebellar ataxia-11 (SCAR11), and a t(1;3) translocation of this gene has been associated with neurodevelopmental



abnormalities. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 4. [provided by RefSeq, Dec 2011],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images