



TELT Polyclonal Antibody

Catalog No	YP-Ab-06273
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	TCAP
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	TELT Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TCAP;CMD1N;TELT;Telethonin (Titin cap protein)
Observed Band	18kD
Calculated Molecular Weight	18kD
Cell Pathway	Cytoplasm, myofibril, sarcomere .
Tissue Specificity	Heart and skeletal muscle.
Function	disease:Defects in TCAP are a cause of cardiomyopathy familial hypertrophic (CMH) [MIM:192600]; also designated FHC or HCM. Familial hypertrophic cardiomyopathy is a hereditary heart disorder characterized by ventricular hypertrophy, which is usually asymmetric and often involves the interventricular septum. The symptoms include dyspnea, syncope, collapse, palpitations, and chest pain. They can be readily provoked by exercise. The disorder has inter- and intrafamilial variability ranging from benign to malignant forms with high risk of cardiac failure and sudden cardiac death.,disease:Defects in TCAP are a cause of limb-girdle muscular dystrophy type 2G (LGMD2G) [MIM:601954]. LGMD2G is an autosomal recessive degenerative myopathy characterized by proximal and distal muscle weakness and atrophy in the limbs, dystrophic changes on muscle biopsy, and absence of telethonin. Cardiac muscle is



Background

titin-cap(TCAP) Homo sapiens Sarcomere assembly is regulated by the muscle protein titin. Titin is a giant elastic protein with kinase activity that extends half the length of a sarcomere. It serves as a scaffold to which myofibrils and other muscle related proteins are attached. This gene encodes a protein found in striated and cardiac muscle that binds to the titin Z1-Z2 domains and is a substrate of titin kinase, interactions thought to be critical to sarcomere assembly. Mutations in this gene are associated with limb-girdle muscular dystrophy type 2G. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night