




TIPRL Polyclonal Antibody

Catalog No	YP-Ab-06282
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	TIPRL
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	TIPRL Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TIP;TIP41;TIP41 like protein;TIPRL;Type 2A interacting protein;TIP41-like protein (Putative MAPK-activating protein PM10) (Type 2A-interacting protein) (TIP)
Observed Band	29kD
Calculated Molecular Weight	29kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Fetal kidney,Lung fibroblast,Muscle,
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Belongs to the TIP41 family.,
Background	TIPRL is an inhibitory regulator of protein phosphatase-2A (PP2A) (see PPP2CA; MIM 176915), PP4 (see PPP4C; MIM 602035), and PP6 (see PPP6C; MIM 612725) (McConnell et al., 2007 [PubMed 17384681]).[supplied by OMIM, Nov 2010],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night