



VGLL2 Polyclonal Antibody

Catalog No	YP-Ab-06289
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	VGLL2 VITO1
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	VGLL2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	VGLL2;VGLL2 VITO1;Transcription cofactor vestigial-like protein 2 (Vgl-2) (Protein VITO1)
Observed Band	34kD
Calculated Molecular Weight	34kD
Cell Pathway	Nucleus .
Tissue Specificity	Skeletal muscle.
Function	function:May act as a specific coactivator for the mammalian TEFs. May play a role in the development of skeletal muscles.,similarity:Belongs to the vestigial family.,subunit:Interacts with TEFs. Binds to TEAD1/TEF1.,tissue specificity:Skeletal muscle.,
Background	This gene encodes a protein with a transcriptional enhancer factor 1 (TEF-1) interaction domain. The encoded protein may act as a co-factor of TEF-1 regulated gene expression during skeletal muscle development. Alternatively spliced transcript variants encoding multiple isoforms have been observed. [provided by RefSeq, Jul 2012],


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from HCT116 cells, primary antibody was diluted at 1:1000, 4° over night