



TAF1L Polyclonal Antibody

Catalog No	YP-Ab-06302
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	TAF1L
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	TAF1L Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TAF(II)210;TAF1L;TAF2A2;TBP associated factor 1 like;TBP associated factor 210 kDa;Transcription initiation factor TFIID subunit 1-like (TAF(II)210) (TBP-associated factor 1-like) (TBP-associated factor 210 kDa) (Transcription initiation factor TFIID 210 kDa subunit)
Observed Band	250-260 kDa
Calculated Molecular Weight	207 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Testis specific, expressed apparently in germ cells.
Function	function:May act as a functional substitute for TAF1/TAFII250 during male meiosis, when sex chromosomes are transcriptionally silenced., similarity:Belongs to the TAF1 family.,similarity:Contains 2 bromo domains., subunit:Can bind directly to TATA-box binding protein (TBP).,tissue specificity:Testis specific, expressed apparently in germ cells.,
Background	This locus is intronless, and apparently arose in the primate lineage from retrotransposition of the transcript from the multi-exon TAF1 locus on the X chromosome. The gene is expressed in male germ cells, and the product has been shown to function interchangeably with the TAF1 product. [provided by

RefSeq, Aug 2015],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images