



TNNT1 Polyclonal Antibody

Catalog No	YP-Ab-06324
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	TNNT1 TNT
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	TNNT1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ANM;STNT;TNNT1;TNT;TNTS;TNNT1 TNT;Troponin T;slow skeletal muscle (TnTs) (Slow skeletal muscle troponin T) (sTnT)
Observed Band	30kD
Calculated Molecular Weight	30kD
Cell Pathway	cytosol,troponin complex,
Tissue Specificity	Skeletal muscle,
Function	disease:Defects in TNNT1 are the cause of nemaline myopathy type 5 (NEM5) [MIM:605355]; also known as Amish nemaline myopathy (ANM) [MIM:605355]. This form of nemaline myopathy (NEM) is common among Old Order Amish with an incidence of approximately 1:500. Affected infants display tremors with hypotonia and mild contractures of the shoulders and hips. Proximal contractures progressively weaken and a pectus carinatum deformity develops before children die of respiratory insufficiency, usually in the second year. A nucleotide replacement in exon 11 causes the protein to be truncated after residue 178.,function:Troponin T is the tropomyosin-binding subunit of troponin, the thin filament regulatory complex which confers calcium-sensitivity to striated muscle actomyosin ATPase activity.,similarity:Belongs to the troponin T family.,



Background

This gene encodes a protein that is a subunit of troponin, which is a regulatory complex located on the thin filament of the sarcomere. This complex regulates striated muscle contraction in response to fluctuations in intracellular calcium concentration. This complex is composed of three subunits: troponin C, which binds calcium, troponin T, which binds tropomyosin, and troponin I, which is an inhibitory subunit. This protein is the slow skeletal troponin T subunit. Mutations in this gene cause nemaline myopathy type 5, also known as Amish nemaline myopathy, a neuromuscular disorder characterized by muscle weakness and rod-shaped, or nemaline, inclusions in skeletal muscle fibers which affects infants, resulting in death due to respiratory insufficiency, usually in the second year. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Ju]


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4° over night