



PTN2 Polyclonal Antibody

Catalog No	YP-Ab-06328
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	PTPN2 PTPT
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	PTN2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	T-cell protein-tyrosine phosphatase;TC PTP;PTPT;PTPN 2;PTN2;PTPN2 PTPT; Tyrosine-protein phosphatase non-receptor type 2 (EC 3.1.3.48) (T-cell protein-tyrosine phosphatase) (TCPTP)
Observed Band	45kD
Calculated Molecular Weight	45kD
Cell Pathway	[Isoform 1]: Endoplasmic reticulum . Endoplasmic reticulum-Golgi intermediate compartment . Targeted to the endoplasmic reticulum by its C-terminal hydrophobic region. .; [Isoform 2]: Nucleus. Cytoplasm. Cell membrane. Predominantly localizes to chromatin (By similarity). Able to shuttle between the nucleus and the cytoplasm and to dephosphorylate plasma membrane receptors (PubMed:9488479). Recruited by activated ITGA1 at the plasma membrane. .
Tissue Specificity	Ubiquitously expressed. Isoform 2 is probably the major isoform. Isoform 1 is expressed in T-cells and in placenta.
Function	catalytic activity:Protein tyrosine phosphate + H(2)O = protein tyrosine + phosphate.,similarity:Belongs to the protein-tyrosine phosphatase family. Non-receptor class 1 subfamily.,similarity:Contains 1 tyrosine-protein phosphatase domain.,subunit:Interacts with FAM82A2.,tissue specificity:PTPA isoform is probably the major PTP expressed in human tissues. PTPB isoform was found



in T-cells and in placenta.,

Background

The protein encoded by this gene is a member of the protein tyrosine phosphatase (PTP) family. Members of the PTP family share a highly conserved catalytic motif, which is essential for the catalytic activity. PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitotic cycle, and oncogenic transformation. Epidermal growth factor receptor and the adaptor protein Shc were reported to be substrates of this PTP, which suggested the roles in growth factor mediated cell signaling. Multiple alternatively spliced transcript variants encoding different isoforms have been found. Two highly related but distinctly processed pseudogenes that localize to chromosomes 1 and 13, respectively, have been reported. [provided by RefSeq, May 2011],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night