



UBP12 Polyclonal Antibody


Catalog No	YP-Ab-06335
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	USP12 UBH1 USP12L1
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	UBP12 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Deubiquitinating enzyme 12;EC:3.4.19.12;UBH1;Ubiquitin carboxyl-terminal hydrolase 12;Ubiquitin hydrolyzing enzyme 1;UBP12;USP12 UBH1 USP12L1; Ubiquitin carboxyl-terminal hydrolase 12 (EC 3.4.19.12) (Deubiquitinating enzyme 12) (Ubiquitin thioesterase 12) (Ubiquitin-hydrolyzing enzyme 1) (Ubiquitin-specific-processing protease 12)
Observed Band	40kD
Calculated Molecular Weight	40kD
Cell Pathway	
Tissue Specificity	Testis,
Function	catalytic activity:Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol., similarity:Belongs to the peptidase C19 family.,
Background	catalytic activity:Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol., similarity:Belongs to the peptidase C19 family.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from SW480 cells,
primary antibody was diluted at 1:1000, 4° over night