



UBP26 Polyclonal Antibody

Catalog No	YP-Ab-06338
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	USP26
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	UBP26 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	UBP26;USP26;Ubiquitin carboxyl-terminal hydrolase 26 (EC 3.4.19.12) (Deubiquitinating enzyme 26) (Ubiquitin thioesterase 26) (Ubiquitin-specific-processing protease 26)
Observed Band	100kD
Calculated Molecular Weight	100kD
Cell Pathway	Nucleus .
Tissue Specificity	PCR rescued clones,Testis,
Function	catalytic activity:Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol., function:Involved in the ubiquitin-dependent proteolytic pathway in conjunction with the 26S proteasome.,similarity:Belongs to the peptidase C19 family.,
Background	This gene encodes a member of the ubiquitin-specific processing (UBP) family of proteases and is a deubiquitinating enzyme (DUB) with His and Cys domains. It is specifically expressed in testis tissue. Mutations in this gene have been associated with Sertoli cell-only syndrome and male infertility. [provided by RefSeq, Jul 2008],

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images