



UBP51 Polyclonal Antibody

Catalog No	YP-Ab-06347
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	USP51
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	UBP51 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	UBP51;USP51;Ubiquitin carboxyl-terminal hydrolase 51 (EC 3.4.19.12) (Deubiquitinating enzyme 51) (Ubiquitin thioesterase 51) (Ubiquitin-specific-processing protease 51)
Observed Band	78kD
Calculated Molecular Weight	78kD
Cell Pathway	Chromosome . Upon DNA damage, chromatin-bound USP51 disassociates from chromatin immediately following DNA damage and reassociates with chromatin following DNA repair. .
Tissue Specificity	Expressed in prostate, brain, lung, aorta and kidney.
Function	catalytic activity:Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol., sequence caution:Contaminating sequence. Potential poly-A sequence., similarity:Belongs to the peptidase C19 family.,similarity:Contains 1 UBP-type zinc finger.,tissue specificity:Expressed in prostate, brain, lung, aorta and kidney.,
Background	catalytic activity:Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol., sequence caution:Contaminating sequence. Potential poly-A sequence., similarity:Belongs to the peptidase C19 family.,similarity:Contains 1 UBP-type

zinc finger.,tissue specificity:Expressed in prostate, brain, lung, aorta and kidney.,

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images