



UBAP2 Polyclonal Antibody

Catalog No	YP-Ab-06355
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	UBAP2 KIAA1491
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	UBAP2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	240890F6;bA176F3.5;KIAA1491;UBAP 2;ubiquitin associated protein 2;UBAP2;UBAP2 KIAA1491;Ubiquitin-associated protein 2 (UBAP-2)
Observed Band	140 kDa
Calculated Molecular Weight	117 kDa
Cell Pathway	nucleus,cytoplasm,cell-cell adherens junction,
Tissue Specificity	Adrenal gland,Brain,Epithelium,Seminom,T-cell,Testis,
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Contains 1 UBA domain.,
Background	The protein encoded by this gene contains a UBA (ubiquitin associated) domain, which is characteristic of proteins that function in the ubiquitination pathway. This gene may show increased expression in the adrenal gland and lymphatic tissues. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2013],
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images