



UB2D4 Polyclonal Antibody

Catalog No	YP-Ab-06358
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	UBE2D4 UBCH5D
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	UB2D4 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	E2(17)KB1;SFT;Stimulator of Fe transport;UBC4/5;UBC4/5 homolog;UBC5A;UBCH5;UBCH5A;UBE2D;UBE2D1;UBE2D1/2/3/4;UBE2D2;UBE2D3;UBE2D4;Ubiquitin carrier protein D1;Ubiquitin protein ligase D1;UB2D4;UBE2D4 UBCH5D;Ubiquitin-conjugating enzyme E2 D4 (EC 6.3.2.19) (HBUCE1) (Ubiquitin carrier protein D4) (Ubiquitin-protein ligase D4)
Observed Band	16kD
Calculated Molecular Weight	16kD
Cell Pathway	
Tissue Specificity	Brain,Ovary,
Function	catalytic activity:ATP + ubiquitin + protein lysine = AMP + diphosphate + protein N-ubiquityllysine.,function:Catalyzes the covalent attachment of ubiquitin to other proteins.,pathway:Protein modification; protein ubiquitination.,similarity:Belongs to the ubiquitin-conjugating enzyme family.,
Background	catalytic activity:ATP + ubiquitin + protein lysine = AMP + diphosphate + protein N-ubiquityllysine.,function:Catalyzes the covalent attachment of ubiquitin to other proteins.,pathway:Protein modification; protein ubiquitination.,similarity:Belongs to the ubiquitin-conjugating enzyme family.,

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

