



UB2Q1 Polyclonal Antibody

Catalog No	YP-Ab-06365
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	UBE2Q1 NICE5 UBE2Q PRO3094
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	UB2Q1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PRO3094;NICE5;NICE 5;GTAP;EC:2.3.2.23;UB2Q1;UBE2Q1 NICE5 UBE2Q PRO3094;Ubiquitin-conjugating enzyme E2 Q1 (EC 6.3.2.19) (Protein NICE-5) (Ubiquitin carrier protein Q1) (Ubiquitin-protein ligase Q1)
Observed Band	46kD
Calculated Molecular Weight	46kD
Cell Pathway	Nucleus . Cell projection, filopodium . Cytoplasm, cytosol .
Tissue Specificity	Widely expressed.
Function	catalytic activity:ATP + ubiquitin + protein lysine = AMP + diphosphate + protein N-ubiquityllysine.,function:Catalyzes the covalent attachment of ubiquitin to other proteins.,pathway:Protein modification; protein ubiquitination.,similarity:Belongs to the ubiquitin-conjugating enzyme family.,tissue specificity:Widely expressed.,
Background	The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes (E1s), ubiquitin-conjugating enzymes (E2s), and ubiquitin-protein ligases (E3s). This gene encodes a member of the E2 ubiquitin-conjugating enzyme family. The

encoded protein is 98% identical to the mouse counterpart. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4°over night