



RPRD2 Polyclonal Antibody

Catalog No	YP-Ab-06374
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	RPRD2 KIAA0460 HSPC099
Immunogen	Synthesized peptide derived from human protein . at AA range: 1370-1450
Specificity	RPRD2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HSPC099;KIAA0460;Regulation of nuclear pre-mRNA domain-containing protein 2;RPRD2;RPRD2 KIAA0460 HSPC099
Observed Band	220 kDa
Calculated Molecular Weight	155 kDa
Cell Pathway	nucleoplasm,DNA-directed RNA polymerase II, holoenzyme,
Tissue Specificity	Brain,Epithelium,Placenta,Salivary gland,Testis,Umb
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Contains 1 CID domain.,
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Contains 1 CID domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

