



ZBT10 Polyclonal Antibody

Catalog No	YP-Ab-06403
Isotype	IgG
Reactivity	Human;Rat
Applications	WB;ELISA
Gene Name	ZBTB10 RINZF RINZFC
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	ZBT10 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RINZF;RINZFC;ZBTB10;Zinc finger protein RIN ZF;ZBT10;ZBTB10 RINZF RINZFC;Zinc finger and BTB domain-containing protein 10 (Zinc finger protein RIN ZF)
Observed Band	65 kDa, 130 kDa
Calculated Molecular Weight	871 aa, 95 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Brain,PCR rescued clones,Prostate,Testis,
Function	function:May be involved in transcriptional regulation.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 2 C2H2-type zinc fingers.,
Background	function:May be involved in transcriptional regulation.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 2 C2H2-type zinc fingers.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images