




CC50A Polyclonal Antibody

Catalog No	YP-Ab-06428
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	TMEM30A C6orf67 CDC50A
Immunogen	Synthesized peptide derived from human protein . at AA range: 110-190
Specificity	CC50A Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CC50A;TMEM30A C6orf67 CDC50A;Cell cycle control protein 50A (Transmembrane protein 30A)
Observed Band	39kD
Calculated Molecular Weight	39kD
Cell Pathway	Membrane ; Multi-pass membrane protein . Cell membrane. Golgi apparatus. Cytoplasmic vesicle, secretory vesicle membrane . Apical cell membrane .
Tissue Specificity	Amygdala,Brain,Liver,Lymph node,Testis,
Function	similarity:Belongs to the CDC50/LEM3 family.,
Background	similarity:Belongs to the CDC50/LEM3 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from DU145 cells,
primary antibody was diluted at 1:1000, 4° over night