



CDK12 Polyclonal Antibody

Catalog No	YP-Ab-06441
Isotype	IgG
Reactivity	Human;Rat
Applications	WB;ELISA
Gene Name	CDK12 CRK7 CRKRS KIAA0904
Immunogen	Synthesized peptide derived from human protein . at AA range: 50-130
Specificity	CDK12 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CrkRS;CRK7;Cell division cycle 2-related protein kinase 7;CDK12;CDC2-related protein kinase 7;CDK12 CRK7 CRKRS KIAA0904;Cyclin-dependent kinase 12 (EC 2.7.11.22) (EC 2.7.11.23) (Cdc2-related kinase;arginine/serine-rich) (CrkRS) (Cell division cycle 2-related protein kinase 7) (CDC2-related protein kinase 7) (Cell divi
Observed Band	180-200 kDa
Calculated Molecular Weight	1490 aa, 164 kDa
Cell Pathway	Nucleus . Nucleus speckle . Colocalized with nuclear speckles throughout interphase. .
Tissue Specificity	Widely expressed.
Function	similarity:Belongs to the protein kinase superfamily.,
Background	similarity:Belongs to the protein kinase superfamily.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

