



CUL7 Polyclonal Antibody

Catalog No	YP-Ab-06457
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	CUL7 KIAA0076
Immunogen	Synthesized peptide derived from human protein . at AA range: 1570-1650
Specificity	CUL7 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CUL 7;CUL7;cullin 7;KIAA0076;CUL7 KIAA0076;Cullin-7 (CUL-7)
Observed Band	185 kDa
Calculated Molecular Weight	1698 aa, 191 kDa
Cell Pathway	Cytoplasm. Cytoplasm, cytoskeleton, microtubule organizing center, centrosome. Cytoplasm, perinuclear region. Golgi apparatus. Colocalizes with FBXW8 at the Golgi apparatus in neurons; localization to Golgi is mediated by OBSL1. During mitosis, localizes to the mitotic apparatus (PubMed:24793695). CCDC8 is required for centrosomal location (PubMed:24793695). .
Tissue Specificity	Highly expressed in fetal kidney and adult skeletal muscle. Also abundant in fetal brain, as well as in adult pancreas, kidney, placenta and heart. Detected in trophoblasts, lymphoblasts, osteoblasts, chondrocytes and skin fibroblasts.
Function	disease:Defects in CUL7 are the cause of 3M syndrome (3MS) [MIM:273750]. 3M syndrome is an autosomal recessive disorder characterized by severe pre- and postnatal growth retardation, facial dysmorphism, large head circumference, and normal intelligence and endocrine function. Skeletal changes include long slender tubular bones and tall vertebral bodies.,function:Component of a probable SCF-like E3 ubiquitin-protein ligase complex, which mediates the ubiquitination and subsequent proteosomal degradation of target proteins.



Probably plays a role in the degradation of proteins involved in endothelial proliferation and/or differentiation (By similarity). Seems not to promote polyubiquitination and proteosomal degradation of TP53. In vitro, complexes of CUL7 with either CUL9 or FBXW8 or TP53 contain E3 ubiquitin-protein ligase activity.,pathway:Protein modification; protein ubiquitination.,si

Background

The protein encoded by this gene is a component of an E3 ubiquitin-protein ligase complex. The encoded protein interacts with TP53, CUL9, and FBXW8 proteins. Defects in this gene are a cause of 3M syndrome type 1 (3M1). Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2009],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images