




MAGBA Polyclonal Antibody

Catalog No	YP-Ab-06475
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	MAGEB10
Immunogen	Synthesized peptide derived from human protein . at AA range: 40-120
Specificity	MAGBA Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MAGBA;MAGEB10;Melanoma-associated antigen B10 (MAGE-B10 antigen)
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	
Tissue Specificity	Blood,Testis,
Function	similarity:Contains 1 MAGE domain.,
Background	MAGE family member B10(MAGEB10) Homo sapiens This gene encodes a member of the B subfamily of the melanoma associated antigen protein family. The encoded protein is specifically expressed in testis and tumor cells. [provided by RefSeq, Apr 2010],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from K562 cells,
primary antibody was diluted at 1:1000, 4°over night