




MAGBI Polyclonal Antibody

Catalog No	YP-Ab-06476
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	MAGEB18
Immunogen	Synthesized peptide derived from human protein . at AA range: 140-220
Specificity	MAGBI Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MAGE B18 antigen;MAGEB18;melanoma antigen family B;18;MAGBI; Melanoma-associated antigen B18 (MAGE-B18 antigen)
Observed Band	37kD
Calculated Molecular Weight	37kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Testis,
Function	similarity:Contains 1 MAGE domain.,
Background	similarity:Contains 1 MAGE domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from SH-SY5Y cells,
primary antibody was diluted at 1:1000, 4° over night