



MAGE2 Polyclonal Antibody

Catalog No	YP-Ab-06480
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	MAGEE2 HCA3
Immunogen	Synthesized peptide derived from human protein . at AA range: 240-320
Specificity	MAGE2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MAGE2;MAGEE2 HCA3;Melanoma-associated antigen E2 (Hepatocellular carcinoma-associated protein 3) (MAGE-E2 antigen)
Observed Band	57kD
Calculated Molecular Weight	57kD
Cell Pathway	
Tissue Specificity	Amygdala,Brain,
Function	similarity:Contains 2 MAGE domains.,
Background	This gene encodes a member of the E subfamily of MAGE (melanoma antigen-encoding gene) gene family. The gene is intronless and the encoded protein has two of the MAGE domains which are characteristic of MAGE family proteins. [provided by RefSeq, Sep 2011],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images