




MAGH1 Polyclonal Antibody

Catalog No	YP-Ab-06482
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	MAGEH1 APR1
Immunogen	Synthesized peptide derived from human protein . at AA range: 10-90
Specificity	MAGH1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Restin;Melanoma-associated antigen H1;MAGE-H1 antigen;APR-1;APR1;MAGH1;MAGEH1 APR1;Melanoma-associated antigen H1 (Apoptosis-related protein 1) (APR-1) (MAGE-H1 antigen) (Restin)
Observed Band	24kD
Calculated Molecular Weight	24kD
Cell Pathway	cytoplasm,
Tissue Specificity	Brain,Lung,Testis,
Function	similarity:Contains 1 MAGE domain.,
Background	This gene belongs to the non-CT (non cancer/testis) subgroup of the melanoma-associated antigen (MAGE) superfamily. The encoded protein is likely associated with apoptosis, cell cycle arrest, growth inhibition or cell differentiation. The protein may be involved in the atRA (all-trans retinoic acid) signaling through the STAT1-alpha (signal transducer and activator of transcription 1-alpha) pathway. [provided by RefSeq, Aug 2013],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4° over night.