



HEN2 Polyclonal Antibody

Catalog No	YP-Ab-06493
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	NHLH2 BHLHA34 HEN2 KIAA0490
Immunogen	Synthesized peptide derived from human protein . at AA range: 30-110
Specificity	HEN2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HEN2;NHLH2 BHLHA34 HEN2 KIAA0490;Helix-loop-helix protein 2 (HEN-2) (Class A basic helix-loop-helix protein 34) (bHLHa34) (Nescient helix loop helix 2) (NSCL-2)
Observed Band	14kD
Calculated Molecular Weight	14kD
Cell Pathway	Nucleus .
Tissue Specificity	Brain,
Function	function:May serve as DNA-binding protein and may be involved in the control of cell-type determination, possibly within the developing nervous system., similarity:Contains 1 basic helix-loop-helix (bHLH) domain.,subunit:Efficient DNA binding requires dimerization with another bHLH protein.,
Background	function:May serve as DNA-binding protein and may be involved in the control of cell-type determination, possibly within the developing nervous system., similarity:Contains 1 basic helix-loop-helix (bHLH) domain.,subunit:Efficient DNA binding requires dimerization with another bHLH protein.,


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night