



DMRTA Polyclonal Antibody

Catalog No	YP-Ab-06505
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	DMRTA1 DMO
Immunogen	Synthesized peptide derived from human protein . at AA range: 170-250
Specificity	DMRTA Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DMRTA;DMRTA1 DMO;Doublesex- and mab-3-related transcription factor A1
Observed Band	55kD
Calculated Molecular Weight	55kD
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in liver, kidney, pancreas, prostate and weakly detected in testis and ovary.
Function	function:May be involved in sexual development.,similarity:Belongs to the DMRT family.,similarity:Contains 1 DM DNA-binding domain.,tissue specificity:Expressed in liver, kidney, pancreas, prostate and weakly detected in testis and ovary.,
Background	function:May be involved in sexual development.,similarity:Belongs to the DMRT family.,similarity:Contains 1 DM DNA-binding domain.,tissue specificity:Expressed in liver, kidney, pancreas, prostate and weakly detected in testis and ovary.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images