



ZNF25 Polyclonal Antibody

Catalog No	YP-Ab-06506
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	ZNF25 KOX19
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	ZNF25 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C2H2 25;CH2H2 25;zinc finger protein 211;Zinc finger protein C2H2 25;ZNF 25;ZNF211;ZNFC25;ZNF25;ZNF25 KOX19;Zinc finger protein 25 (Zinc finger protein KOX19)
Observed Band	60 kDa
Calculated Molecular Weight	576 aa, 66 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Hippocampus,Lymphoid,Teratocarcinoma,Testis,
Function	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 12 C2H2-type zinc fingers.,
Background	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 12 C2H2-type zinc fingers.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images